

— ABLYBIO, Help Your Research



# M-CSF (YD18249) Rabbit mAb

货号: **AYD14938**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	60kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,Jurkat,K-562
细胞定位	Cell membrane, Secreted, extracellular space
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants.
基因ID	1435
基因名	CSF1
Swiss	P09603 ( <a href="https://www.uniprot.org/uniprotkb/P09603/entry">https://www.uniprot.org/uniprotkb/P09603/entry</a> )
别名	M-CSF (YD18249),M-CSF (YD18249) Rabbit mAb,CSF1,Lanimostim,Proteoglycan macrophage colony-stimulating factor

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)