

— ABLYBIO, Help Your Research



Osteopontin (YD15880) Rabbit mAb

货号: **AYD14968**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney
细胞定位	Secreted
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is involved in the attachment of osteoclasts to the mineralized bone matrix. The encoded protein is secreted and binds hydroxyapatite with high affinity. The osteoclast vitronectin receptor is found in the cell membrane and may be involved in the binding to this protein. This protein is also a cytokine that upregulates expression of interferon-gamma and interleukin-12. Several transcript variants encoding different isoforms have been found for this gene.
基因ID	6696
基因名	SPP1
Swiss	P10451 (https://www.uniprot.org/uniprotkb/P10451/entry)
别名	Osteopontin (YD15880),Osteopontin (YD15880) Rabbit mAb,SPP1,Bone sialoprotein 1,Nephropontin,Secreted phosphoprotein 1,Urinary stone protein,Uropontin,BNSP,OPN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)