

— ABLYBIO, Help Your Research



PDX1 (YD15816) Rabbit mAb

货号: **AYD14986**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse intestine
细胞定位	Nucleus, Cytoplasm, cytosol
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a transcriptional activator of several genes, including insulin, somatostatin, glucokinase, islet amyloid polypeptide, and glucose transporter type 2. The encoded nuclear protein is involved in the early development of the pancreas and plays a major role in glucose-dependent regulation of insulin gene expression. Defects in this gene are a cause of pancreatic agenesis, which can lead to early-onset insulin-dependent diabetes mellitus (NIDDM), as well as maturity onset diabetes of the young type 4 (MODY4).
基因ID	3651
基因名	PDX1
Swiss	P52945 (https://www.uniprot.org/uniprotkb/P52945/entry)
别名	PDX1 (YD15816),PDX1 (YD15816) Rabbit mAb,PDX1,Glucose-sensitive factor,Insulin promoter factor 1,Insulin upstream factor 1,Islet/duodenum homeobox-1,Somatostatin-transactivating factor 1,IPF1,STF1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)