

— ABLYBIO, Help Your Research



c-Fos (YD15746) Rabbit mAb

货号: **AYD14998**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	41kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus, Endoplasmic reticulum, Cytoplasm, cytosol
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. In some cases, expression of the FOS gene has also been associated with apoptotic cell death.
基因ID	2353
基因名	FOS
Swiss	P01100 (https://www.uniprot.org/uniprotkb/P01100/entry)
别名	c-Fos (YD15746),c-Fos (YD15746) Rabbit mAb,FOS,Cellular oncogene fos,Fos proto-oncogene,AP-1 transcription factor subunit,G0/G1 switch regulatory protein 7,Proto-oncogene c-Fos,Transcription factor AP-1 subunit c-Fos,G0S7

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)