

Myeloperoxidase (YD15738) Rabbit mAb

货号: **AYD15000**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60
细胞定位	Lysosome
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	Myeloperoxidase (MPO) is a heme protein synthesized during myeloid differentiation that constitutes the major component of neutrophil azurophilic granules. Produced as a single chain precursor, myeloperoxidase is subsequently cleaved into a light and heavy chain. The mature myeloperoxidase is a tetramer composed of 2 light chains and 2 heavy chains. This enzyme produces hypohalous acids central to the microbicidal activity of neutrophils.
基因ID	4353
基因名	MPO
Swiss	P05164

别名	Myeloperoxidase (YD15738)
----	---------------------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn