

— ABLYBIO, Help Your Research



CD69 (YD15720) Rabbit mAb

货号: **AYD15006**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-937,Mouse lung,Mouse spleen,Rat lung,Rat spleen
细胞定位	Cell membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the calcium dependent lectin superfamily of type II transmembrane receptors. Expression of the encoded protein is induced upon activation of T lymphocytes, and may play a role in proliferation. Furthermore, the protein may act to transmit signals in natural killer cells and platelets.
------	---

基因ID	969
基因名	CD69
Swiss	Q07108 (https://www.uniprot.org/uniprotkb/Q07108/entry)
别名	CD69 (YD15720),CD69 (YD15720) Rabbit mAb,CD69,Activation inducer molecule,BL-AC/P26,C-type lectin domain family 2 member C,EA1,Early T-cell activation antigen p60,GP32/28,Leukocyte surface antigen Le u-23,MLR-3,CLEC2C

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)