

— ABLYBIO, Help Your Research



# MyoD1 (YD15696) Rabbit mAb

货号: **AYD15013**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	35kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse skeletal muscle
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a nuclear protein that belongs to the basic helix-loop-helix family of transcription factors and the myogenic factors subfamily. It regulates muscle cell differentiation by inducing cell cycle arrest, a prerequisite for myogenic initiation. The protein is also involved in muscle regeneration. It activates its own transcription which may stabilize commitment to myogenesis.
------	--

基因ID	4654
基因名	MYOD1
Swiss	P15172 ( <a href="https://www.uniprot.org/uniprotkb/P15172/entry">https://www.uniprot.org/uniprotkb/P15172/entry</a> )
别名	MyoD1 (YD15696),MyoD1 (YD15696) Rabbit mAb,MYOD1,Class C basic helix-loop-helix protein 1,Myogenic factor 3,BHLHC1,MYF3,MYOD

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)