

— ABLYBIO, Help Your Research



# N Cadherin (YD15541) Rabbit mAb

货号: **AYD15040**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	100kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Cell membrane, sarcolemma, Cell junction, Cell surface, desmosome, adherens junction, Cell projection, axon
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a classical cadherin and member of the cadherin superfamily. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein is proteolytically processed to generate a calcium-dependent cell adhesion molecule and glycoprotein. This protein plays a role in the establishment of left-right asymmetry, development of the nervous system and the formation of cartilage and bone.
基因ID	1000
基因名	CDH2
Swiss	P19022 ( <a href="https://www.uniprot.org/uniprotkb/P19022/entry">https://www.uniprot.org/uniprotkb/P19022/entry</a> )
别名	N Cadherin (YD15541),N Cadherin (YD15541) Rabbit mAb,CDH2,CDw325,Neural cadherin,CDHN,NCAD

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)