

# TIMP1 (YD15448) Rabbit mAb

货号: **AYD15052**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	

## 抗原信息

抗原信息	
------	--

## 靶点信息


研究背景	Metalloproteinase inhibitor that functions by forming one to one complexes with target metalloproteinase s, such as collagenases, and irreversibly inactivates them by binding to their catalytic zinc cofactor. Acts on MMP1, MMP2, MMP3, MMP7, MMP8, MMP9, MMP10, MMP11, MMP12, MMP13 and MMP16. Does not act on MMP14 (By similarity). Also functions as a growth factor that regulates cell differentiation, migration and cell death and activates cellular signaling cascades via CD63 and ITGB1. Plays a role in integrin signaling
基因ID	3134
基因名	Timp1
Swiss	P12032

---

别名	TIMP1 (YD15448)
----	-----------------

### 产品验证

---



### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)