

— ABLYBIO, Help Your Research



# HMGA1 (YD15351) Rabbit mAb

货号: **AYD15083**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Nucleus, Chromosome
纯化	亲和纯化




## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a chromatin-associated protein involved in the regulation of gene transcription, integration of retroviruses into chromosomes, and the metastatic progression of cancer cells. The encoded protein preferentially binds to the minor groove of AT-rich regions in double-stranded DNA. Multiple transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene have been identified on multiple chromosomes.
基因ID	3159
基因名	HMGA1
Swiss	P17096 ( <a href="https://www.uniprot.org/uniprotkb/P17096/entry">https://www.uniprot.org/uniprotkb/P17096/entry</a> )
别名	HMGA1 (YD15351),HMGA1 (YD15351) Rabbit mAb,HMGA1,High mobility group AT-hook protein 1,High mobility group protein R,HMG1Y

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)