

— ABLYBIO, Help Your Research



CD8 alpha (YD15321) Rabbit mAb

货号: **AYD15089**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC IP
推荐浓度	
理论分子量	27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen,Mouse thymus,Mouse brain(Negative control)
细胞定位	Cell membrane
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Enables identical protein binding activity. Acts upstream of or within several processes, including cytotoxic T cell differentiation; defense response to virus; and positive regulation of calcium-mediated signaling. Located in external side of plasma membrane. Is expressed in several structures, including endocrine gland; exocrine gland; genitourinary system; gut; and hemolymphoid system gland. Orthologous to human CD8A (CD8a molecule).
基因ID	12525
基因名	Cd8a
Swiss	P01731 (https://www.uniprot.org/uniprotkb/P01731/entry)
别名	CD8 alpha (YD15321),CD8 alpha (YD15321) Rabbit mAb,Cd8a,T-cell surface glycoprotein Lyt-2,Lyt-2,Lyt2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)