

— ABLYBIO, Help Your Research



# IL-1 beta (YD15299) Rabbit mAb

货号: **AYD15096**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	Cytoplasm, cytosol, Secreted, Lysosome, extracellular exosome
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1 (CASP1/ICE). This cytokine is an important mediator of the inflammatory response, and is involved in a variety of cellular activities, including cell proliferation, differentiation, and apoptosis. The induction of cyclooxygenase-2 (PTGS2/COX2) by this cytokine in the central nervous system (CNS) is found to contribute to inflammatory pain hypersensitivity. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2.
基因ID	3553
基因名	IL1B
Swiss	P01584 ( <a href="https://www.uniprot.org/uniprotkb/P01584/entry">https://www.uniprot.org/uniprotkb/P01584/entry</a> )
别名	IL-1 beta (YD15299),IL-1 beta (YD15299) Rabbit mAb,IL1B,Catabolin,IL1F2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)