

— ABLYBIO, Help Your Research



# Gephyrin (YD15193) Rabbit mAb

货号: **AYD15115**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	80kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse brain,Mouse liver,Rat brain
细胞定位	Postsynaptic cell membrane, Cell membrane, Cytoplasm, cytosol, cytoskeleton, Cell projection, dendrite, Postsynaptic density
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a neuronal assembly protein that anchors inhibitory neurotransmitter receptors to the postsynaptic cytoskeleton via high affinity binding to a receptor subunit domain and tubulin dimers. In non-neuronal tissues, the encoded protein is also required for molybdenum cofactor biosynthesis. Mutations in this gene may be associated with the neurological condition hyperplexia and also lead to molybdenum cofactor deficiency. Numerous alternatively spliced transcript variants encoding different isoforms have been described; however, the full-length nature of all transcript variants is not currently known.
基因ID	10243
基因名	GPHN
Swiss	Q9NQX3 ( <a href="https://www.uniprot.org/uniprotkb/Q9NQX3/entry">https://www.uniprot.org/uniprotkb/Q9NQX3/entry</a> )
别名	Gephyrin (YD15193), Gephyrin (YD15193) Rabbit mAb, GPHN, Domain G, Domain E, Mo insertase, GPH, KIAA1385

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)