

— ABLYBIO, Help Your Research



# Neuroigin 4 (YD15128) Rabbit mAb

货号: **AYD15131**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	92kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,Mouse brain,Mouse thymus,Rat brain
细胞定位	Cell membrane, Postsynaptic density membrane, Cell projection, dendritic spine, dendrite, Synapse
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the type-B carboxylesterase/lipase protein family. The encoded protein belongs to a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. The encoded protein interacts with discs large homolog 4 (DLG4). Mutations in this gene have been associated with autism and Asperger syndrome. Alternative splicing results in multiple transcript variants.
基因ID	57502
基因名	NLGN4X
Swiss	Q8N0W4 ( <a href="https://www.uniprot.org/uniprotkb/Q8N0W4/entry">https://www.uniprot.org/uniprotkb/Q8N0W4/entry</a> )
别名	Neurologin 4 (YD15128),Neurologin 4 (YD15128) Rabbit mAb,NLGN4X,HNLX,KIAA1260,NLGN4

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)