

— ABLYBIO, Help Your Research



# BCRP/ABCG2 (YD15077) Rabbit mAb

货号: **AYD15150**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Cell membrane, Apical cell membrane, Mitochondrion membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The membrane-associated protein encoded by this gene is included in the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. Alternatively referred to as a breast cancer resistance protein, this protein functions as a xenobiotic transporter which may play a major role in multi-drug resistance. It likely serves as a cellular defense mechanism in response to mitoxantrone and anthracycline exposure. Significant expression of this protein has been observed in the placenta, which may suggest a potential role for this molecule in placenta tissue. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	9429
基因名	ABCG2
Swiss	Q9UNQ0 ( <a href="https://www.uniprot.org/uniprotkb/Q9UNQ0/entry">https://www.uniprot.org/uniprotkb/Q9UNQ0/entry</a> )
别名	BCRP/ABCG2 (YD15077),BCRP/ABCG2 (YD15077) Rabbit mAb,ABCG2,ATP-binding cassette sub-family G member 2,Breast cancer resistance protein,CDw338,Mitoxantrone resistance-associated protein,Placenta-specific ATP-binding cassette transporter,Urate exporter,ABCP,BCRP,BCRP1,MXR

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)