

— ABLYBIO, Help Your Research



FGF21 (YD15069) Rabbit mAb

货号: **AYD15154**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse pancreas
细胞定位	Secreted
纯化	亲和纯化

抗原信息



抗原信息	
------	--

靶点信息

研究背景	Theis gene encodes a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes. This protein is a secreted endocrine factor that functions as a major metabolic regulator. The encoded protein stimulates the uptake of glucose in adipose tissue.
------	---

基因ID	26291
基因名	FGF21
Swiss	Q9NSA1 (https://www.uniprot.org/uniprotkb/Q9NSA1/entry)
别名	FGF21 (YD15069),FGF21 (YD15069) Rabbit mAb,FGF21

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)