

— ABLYBIO, Help Your Research



Syk (YD14958) Rabbit mAb

货号: **AYD15187**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P ICC/IF |
| 推荐浓度 | |
| 理论分子量 | 72kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Jurkat,Mouse spleen,Rat liver |
| 细胞定位 | Cell membrane, Cytoplasm, cytosol |
| 纯化 | 亲和纯化 |


抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 6850 |
| 基因名 | SYK |
| Swiss | P43405 (https://www.uniprot.org/uniprotkb/P43405/entry) |
| 别名 | Syk (YD14958), Syk (YD14958) Rabbit mAb, SYK, Spleen tyrosine kinase, p72-Syk |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)