

Bok (YD14877) Rabbit mAb

货号: **AYD15214**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat testis
细胞定位	Mitochondrion membrane, Endoplasmic reticulum membrane, Mitochondrion inner membrane, Cytoplasm , Nucleus, Mitochondrion, Endoplasmic reticulum, Mitochondrion outer membrane, Early endosome membrane, Recycling endosome membrane, Nucleus outer membrane, Golgi apparatus, cis-Golgi network membrane, trans-Golgi network membrane, Membrane
纯化	

抗原信息

抗原信息	
------	--


靶点信息

研究背景	The protein encoded by this gene belongs to the BCL2 family, members of which form homo- or heterodimers, and act as anti- or proapoptotic regulators that are involved in a wide variety of cellular processes. Studies in rat show that this protein has restricted expression in reproductive tissues, interacts strongly with some antiapoptotic BCL2 proteins, not at all with proapoptotic BCL2 proteins, and induces apoptosis in transfected cells. Thus, this protein represents a proapoptotic member of the BCL2 family.
基因ID	666

基因名	BOK
Swiss	Q9UMX3
别名	Bok (YD14877)

产品验证

Mouse small intestine



实验步骤

访问官网浏览详情: www.ablybio.cn