

— ABLYBIO, Help Your Research



TREM2 (YD14876) Rabbit mAb

货号: **AYD15215**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cell membrane, Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a membrane protein that forms a receptor signaling complex with the TYRO protein tyrosine kinase binding protein. The encoded protein functions in immune response and may be involved in chronic inflammation by triggering the production of constitutive inflammatory cytokines. Defects in this gene are a cause of polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (PLO SL). Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	54209
基因名	TREM2
Swiss	Q9NZC2 (https://www.uniprot.org/uniprotkb/Q9NZC2/entry)
别名	TREM2 (YD14876),TREM2 (YD14876) Rabbit mAb,TREM2,Triggering receptor expressed on monocytes 2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)