

— ABLYBIO, Help Your Research



Protamine 2 (YD14799) Rabbit mAb

货号: **AYD15227**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	13kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus, Chromosome
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Protamines substitute for histones in the chromatin of sperm during the haploid phase of spermatogenesis. They compact sperm DNA into a highly condensed, stable and inactive complex
基因ID	5620

基因名	PRM2
Swiss	P04554 (https://www.uniprot.org/uniprotkb/P04554/entry)
别名	Protamine 2 (YD14799),Protamine 2 (YD14799) Rabbit mAb,PRM2,Sperm histone P2,Sperm protamine P2,Sperm protamine P4,P2'',Sperm protamine P3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)