

— ABLYBIO, Help Your Research



# Intra Acrosomal Protein (YD14778) Rabbit mAb

货号: **AYD15230**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Cytoplasmic vesicle, secretory vesicle, acrosome
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a testis-specific, differentiation antigen, acrosomal vesicle protein 1, that arises within the acrosomal vesicle during spermatogenesis, and is associated with the acrosomal membranes and matrix of mature sperm. The acrosomal vesicle protein 1 may be involved in sperm-zona binding or penetration. Alternatively spliced transcript variants have been described.
------	---

基因ID	56
基因名	ACRV1
Swiss	P26436 ( <a href="https://www.uniprot.org/uniprotkb/P26436/entry">https://www.uniprot.org/uniprotkb/P26436/entry</a> )
别名	Intra Acrosomal Protein (YD14778),Intra Acrosomal Protein (YD14778) Rabbit mAb,ACRV1,Acrosomal vesicle protein 1

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)