

— ABLYBIO, Help Your Research



PHF1 (YD14732) Rabbit mAb

货号: **AYD15245**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	62kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,Mouse lung
细胞定位	Nucleus, Cytoplasm, cytoskeleton, microtubule organizing center, centrosome
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a Polycomb group protein. The protein is a component of a histone H3 lysine-27 (H3K27)-specific methyltransferase complex, and functions in transcriptional repression of homeotic genes. The protein is also recruited to double-strand breaks, and reduced protein levels results in X-ray sensitivity and increased homologous recombination. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	5252
基因名	PHF1
Swiss	O43189 (https://www.uniprot.org/uniprotkb/O43189/entry)
别名	PHF1 (YD14732),PHF1 (YD14732) Rabbit mAb,PHF1,Polycomb-like protein 1,PCL1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)