

AMCase (YD14559) Rabbit mAb

货号: **AYD15272**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HepG2,Mouse stomach,Rat lung,Rat stomach
细胞定位	Secreted, Cytoplasm
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene degrades chitin, which is found in the cell wall of most fungi as well as in arthropods and some nematodes. The encoded protein can also stimulate interleukin 13 expression, and variations in this gene can lead to asthma susceptibility. Several transcript variants encoding a few different isoforms have been found for this gene.
基因ID	27159
基因名	CHIA
Swiss	Q9BZP6

别名	AMCase (YD14559)
----	------------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn