

— ABLYBIO, Help Your Research



BAFF (YD14537) Rabbit mAb

货号: **AYD15276**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-937,Mouse spleen,Rat lung
细胞定位	Cell membrane, Secreted
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B/TACI, TNFRSF17/BCMA, and TNFRSF13C/BAFFR. This cytokine is expressed in B cell lineage cells, and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. Alternatively spliced transcript variants encoding distinct isoforms have been identified.
基因ID	10673
基因名	TNFSF13B
Swiss	Q9Y275 (https://www.uniprot.org/uniprotkb/Q9Y275/entry)
别名	BAFF (YD14537),BAFF (YD14537) Rabbit mAb,TNFSF13B,B lymphocyte stimulator,B-cell-activating factor, BAFF,Dendritic cell-derived TNF-like molecule,TNF- and APOL-related leukocyte expressed ligand 1,BLYS,T ALL1,TNFSF20,ZTNF4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)