

— ABLYBIO, Help Your Research



# ATF3 (YD14528) Rabbit mAb

货号: **AYD15278**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF IP
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,RAW264.7
细胞定位	Nucleus
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the mammalian activation transcription factor/cAMP responsive element-binding (CREB) protein family of transcription factors. This gene is induced by a variety of signals, including many of those encountered by cancer cells, and is involved in the complex process of cellular stress response. Multiple transcript variants encoding different isoforms have been found for this gene. It is possible that alternative splicing of this gene may be physiologically important in the regulation of target genes.
基因ID	467
基因名	ATF3
Swiss	P18847 ( <a href="https://www.uniprot.org/uniprotkb/P18847/entry">https://www.uniprot.org/uniprotkb/P18847/entry</a> )
别名	ATF3 (YD14528),ATF3 (YD14528) Rabbit mAb,ATF3,Activating transcription factor 3

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)