

— ABLYBIO, Help Your Research



CTAG1B (YD14380) Rabbit mAb

货号: **AYD15320**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	18kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080, Jurkat, HeLa
细胞定位	Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is an antigen that is overexpressed in many cancers but that is also expressed in normal testis. This gene is found in a duplicated region of the X-chromosome and therefore has a neighboring gene of identical sequence.
------	---

基因ID	1485
基因名	CTAG1A, CTAG1B
Swiss	P78358 (https://www.uniprot.org/uniprotkb/P78358/entry)
别名	CTAG1B (YD14380),CTAG1B (YD14380) Rabbit mAb,CTAG1A,CTAG1B,Autoimmunogenic cancer/testis antigen NY-ESO-1,Cancer/testis antigen 6.1,L antigen family member 2,CTAG,CTAG1,ESO1,LAGE2,LAGE2A,LAGE2B

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)