

— ABLYBIO, Help Your Research



CD73 (YD14366) Rabbit mAb

货号: **AYD15327**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	63kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A375,Mouse kidney,Rat eye
细胞定位	Cell membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4907
基因名	NT5E
Swiss	P21589 (https://www.uniprot.org/uniprotkb/P21589/entry)
别名	CD73 (YD14366),CD73 (YD14366) Rabbit mAb,NT5E,5'-deoxynucleotidase,Ecto-5'-nucleotidase,IMP-specific 5'-nucleotidase,Thymidylate 5'-phosphatase,NT5,NTE

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)