

— ABLYBIO, Help Your Research



gamma Actin (YD14302) Rabbit mAb

货号: **AYD15346**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm, cytoskeleton
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells. May play a role in the repair of noise-induced stereocilia gaps thereby maintains hearing sensitivity following loud noise damage (By similarity)
------	--

基因ID	71
基因名	ACTG1
Swiss	P63261 (https://www.uniprot.org/uniprotkb/P63261/entry)
别名	gamma Actin (YD14302),gamma Actin (YD14302) Rabbit mAb,ACTG1,Gamma-actin,ACTG

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)