

— ABLYBIO, Help Your Research



# Carbonic Anhydrase II (YD14201) Rabbit mAb

货号: **AYD15388**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,293T,Mouse kidney,Mouse spleen,Mouse lung
细胞定位	Cytoplasm, Cell membrane
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is one of several isozymes of carbonic anhydrase, which catalyzes reversible hydration of carbon dioxide. Defects in this enzyme are associated with osteopetrosis and renal tubular acidosis. Two transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	760
基因名	CA2
Swiss	P00918 ( <a href="https://www.uniprot.org/uniprotkb/P00918/entry">https://www.uniprot.org/uniprotkb/P00918/entry</a> )
别名	Carbonic Anhydrase II (YD14201),Carbonic Anhydrase II (YD14201) Rabbit mAb,CA2,Carbonate dehydratase II,Carbonic anhydrase C,Carbonic anhydrase II,Cyanamide hydratase CA2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)