

— ABLYBIO, Help Your Research



Albumin (YD14128) Rabbit mAb

货号: **AYD15418**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,HeLa(Negative control),Mouse liver,Mouse lung,Rat liver
细胞定位	Secreted
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	Albumin is a soluble, monomeric protein which comprises about one-half of the blood serum protein. Albumin functions primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Albumin is a globular unglycosylated serum protein of molecular weight 65,000. Albumin is synthesized in the liver as prealbumin which has an N-terminal peptide that is removed before the nascent protein is released from the rough endoplasmic reticulum. The product, proalbumin, is in turn cleaved in the Golgi vesicles to produce the secreted albumin.
基因ID	213
基因名	ALB
Swiss	P02768 (https://www.uniprot.org/uniprotkb/P02768/entry)
别名	Albumin (YD14128),Albumin (YD14128) Rabbit mAb,ALB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)