

Steroidogenic Factor 1 (YD14119) Rabbit mAb

货号: **AYD15421**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,Mouse testis
细胞定位	Nucleus
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a transcriptional activator involved in sex determination. The encoded protein binds DNA as a monomer. Defects in this gene are a cause of XY sex reversal with or without adrenal failure as well as adrenocortical insufficiency without ovarian defect.
基因ID	2516
基因名	NR5A1
Swiss	Q13285
别名	Steroidogenic Factor 1 (YD14119)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn