

— ABLYBIO, Help Your Research



Uteroglobin (YD14062) Rabbit mAb

货号: **AYD15441**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	10kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Secreted
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.
------	--

基因ID	7356
基因名	SCGB1A1
Swiss	P11684 (https://www.uniprot.org/uniprotkb/P11684/entry)
别名	Uteroglobin (YD14062),Uteroglobin (YD14062) Rabbit mAb,SCGB1A1,Club cell phospholipid-binding protein,Club cells 10 kDa secretory protein,Secretoglobin family 1A member 1,Urinary protein 1,CC10,CCSP,UG B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)