

— ABLYBIO, Help Your Research



# LFA3 (YD14052) Rabbit mAb

货号: **AYD15444**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P
推荐浓度	
理论分子量	28kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,293T
细胞定位	Cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the immunoglobulin superfamily. The encoded protein is a ligand of the T lymphocyte CD2 protein, and functions in adhesion and activation of T lymphocytes. The protein is localized to the plasma membrane. Alternatively spliced transcript variants have been described.
------	---

基因ID	965
基因名	CD58
Swiss	P19256 ( <a href="https://www.uniprot.org/uniprotkb/P19256/entry">https://www.uniprot.org/uniprotkb/P19256/entry</a> )
别名	LFA3 (YD14052),LFA3 (YD14052) Rabbit mAb,CD58,Surface glycoprotein LFA-3,LFA3

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)