

— ABLYBIO, Help Your Research



# IREB2 (YD14047) Rabbit mAb

货号: **AYD15447**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	105kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Rat liver
细胞定位	Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	RNA-binding protein that binds to iron-responsive elements (IRES, which are stem-loop structures found in the 5'-UTR of ferritin, and delta aminolevulinic acid synthase mRNAs, and in the 3'-UTR of transferrin receptor mRNA. Binding to the IRE element in ferritin results in the repression of its mRNA translation. Binding of the protein to the transferrin receptor mRNA inhibits the degradation of this otherwise rapidly degraded mRNA.
基因ID	3658
基因名	IREB2
Swiss	P48200 ( <a href="https://www.uniprot.org/uniprotkb/P48200/entry">https://www.uniprot.org/uniprotkb/P48200/entry</a> )
别名	IREB2 (YD14047),IREB2 (YD14047) Rabbit mAb,IREB2,Iron regulatory protein 2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)