

— ABLYBIO, Help Your Research



# Versican (YD14034) Rabbit mAb

货号: **AYD15452**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	373kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse brain
细胞定位	Secreted, extracellular space, extracellular matrix, Cell projection, cilium, photoreceptor outer segment, interphotoreceptor matrix
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene is a member of the aggrecan/versican proteoglycan family. The protein encoded is a large chondroitin sulfate proteoglycan and is a major component of the extracellular matrix. This protein is involved in cell adhesion, proliferation, proliferation, migration and angiogenesis and plays a central role in tissue morphogenesis and maintenance. Mutations in this gene are the cause of Wagner syndrome type 1. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2009]
基因ID	1462
基因名	VCAN
Swiss	P13611 ( <a href="https://www.uniprot.org/uniprotkb/P13611/entry">https://www.uniprot.org/uniprotkb/P13611/entry</a> )
别名	Versican (YD14034),Versican (YD14034) Rabbit mAb,VCAN,Chondroitin sulfate proteoglycan core protein 2,Glial hyaluronate-binding protein,Large fibroblast proteoglycan,PG-M,CSPG2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)