

— ABLYBIO, Help Your Research



Hsp105 (YD14032) Rabbit mAb

货号: **AYD15453**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	97kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,BT-474,HepG2,293T,U-251MG,Mouse brain,Mouse kidney,Rat brain
细胞定位	Cytoplasm
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the heat shock protein 70 family of proteins. The encoded protein functions as a nucleotide exchange factor for the molecular chaperone heat shock cognate 71 kDa protein (Hsc70). In addition, this protein plays a distinct but related role as a holdase that inhibits the aggregation of misfolded proteins, including the cystic fibrosis transmembrane conductance regulator (CFTR) protein. Elevated expression of this protein has been observed in numerous human cancers.
基因ID	10808
基因名	HSPH1
Swiss	Q92598 (https://www.uniprot.org/uniprotkb/Q92598/entry)
别名	Hsp105 (YD14032),Hsp105 (YD14032) Rabbit mAb,HSPH1,Antigen NY-CO-25,Heat shock 110 kDa protein, Heat shock protein family H member 1,HSP105,HSP110,KIAA0201

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)