

— ABLYBIO, Help Your Research



MIF (YD13814) Rabbit mAb

货号: **AYD15493**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	12kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HL-60,Mouse liver,Mouse brain,Mouse kidney,Rat liver,Rat brain
细胞定位	Secreted, Cytoplasm
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a lymphokine involved in cell-mediated immunity, immunoregulation, and inflammation. It plays a role in the regulation of macrophage function in host defense through the suppression of anti-inflammatory effects of glucocorticoids. This lymphokine and the JAB1 protein form a complex in the cytosol near the peripheral plasma membrane, which may indicate an additional role in integrin signaling pathways.
基因ID	4282
基因名	MIF
Swiss	P14174 (https://www.uniprot.org/uniprotkb/P14174/entry)
别名	MIF (YD13814),MIF (YD13814) Rabbit mAb,MIF,Glycosylation-inhibiting factor,L-dopachrome isomerase,L-dopachrome tautomerase,Phenylpyruvate tautomerase,GLIF,MMIF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)