

APG5L/ATG5 (YD13696) Rabbit mAb

货号: **AYD15511**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC
推荐浓度	
理论分子量	32kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,NIH/3T3,K-562,C6,SH-SY5Y
细胞定位	Cytoplasm, Preautophagosomal structure membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene, in combination with autophagy protein 12, functions as an E1-like activating enzyme in a ubiquitin-like conjugating system. The encoded protein is involved in several cellular processes, including autophagic vesicle formation, mitochondrial quality control after oxidative damage, negative regulation of the innate antiviral immune response, lymphocyte development and proliferation, MHC II antigen presentation, adipocyte differentiation, and apoptosis. Several transcript variants encoding different protein isoforms have been found for this gene.
基因ID	9474
基因名	ATG5

Swiss	Q9H1Y0
别名	APG5L/ATG5 (YD13696)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn