

DKK1 (YD13687) Rabbit mAb

货号: **AYD15512**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,K-562,Mouse brain,Rat brain,Rat liver
细胞定位	Secreted
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.
基因ID	22943
基因名	DKK1
Swiss	O94907

别名	DKK1 (YD13687)
----	----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn