

— ABLYBIO, Help Your Research



# Bid (YD13684) Rabbit mAb

货号: **AYD15514**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Rat thymus
细胞定位	Cytoplasm, Mitochondrion membrane, Mitochondrion outer membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a death agonist that heterodimerizes with either agonist BAX or antagonist BCL2. The encoded protein is a member of the BCL-2 family of cell death regulators. It is a mediator of mitochondria l damage induced by caspase-8 (CASP8); CASP8 cleaves this encoded protein, and the COOH-terminal part translocates to mitochondria where it triggers cytochrome c release. Multiple alternatively spliced transcript variants have been found, but the full-length nature of some variants has not been defined.
基因ID	637
基因名	BID
Swiss	P55957 ( <a href="https://www.uniprot.org/uniprotkb/P55957/entry">https://www.uniprot.org/uniprotkb/P55957/entry</a> )
别名	Bid (YD13684),Bid (YD13684) Rabbit mAb,BID,p22 BID,Cytochrome c Releasing Factor,Truncated Bid,p13 BID,p11 BID

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)