

— ABLYBIO, Help Your Research



# RAP1A/RAP1B (YD13521) Rabbit mAb

货号: **AYD15548**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG,HL-60
细胞定位	Cell membrane, Cytoplasm, perinuclear region, Cell junction, Early endosome
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	<p>This gene encodes a member of the Ras family of small GTPases. The encoded protein undergoes a change in conformational state and activity, depending on whether it is bound to GTP or GDP. This protein is activated by several types of guanine nucleotide exchange factors (GEFs), and inactivated by two groups of GTPase-activating proteins (GAPs). The activation status of the encoded protein is therefore affected by the balance of intracellular levels of GEFs and GAPs. The encoded protein regulates signaling pathways that affect cell proliferation and adhesion, and may play a role in tumor malignancy. Pseudogenes of this gene have been defined on chromosomes 14 and 17. Alternative splicing results in multiple transcript variants.</p>
基因ID	5906
基因名	RAP1A
Swiss	P62834 ( <a href="https://www.uniprot.org/uniprotkb/P62834/entry">https://www.uniprot.org/uniprotkb/P62834/entry</a> )
别名	RAP1A/RAP1B (YD13521),RAP1A/RAP1B (YD13521) Rabbit mAb,RAP1A,C21KG,G-22K,GTP-binding protein smg p21A,Ras-related protein Krev-1,KREV1

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)