

— ABLYBIO, Help Your Research



Apolipoprotein H (YD13515) Rabbit mAb

货号: **AYD15550**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	38kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1
细胞定位	Secreted
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	Apolipoprotein H has been implicated in a variety of physiologic pathways including lipoprotein metabolism, coagulation, and the production of antiphospholipid autoantibodies. APOH may be a required cofactor for anionic phospholipid binding by the antiphospholipid autoantibodies found in sera of many patients with lupus and primary antiphospholipid syndrome, but it does not seem to be required for the reactivity of antiphospholipid autoantibodies associated with infections.
基因ID	350
基因名	APOH
Swiss	P02749 (https://www.uniprot.org/uniprotkb/P02749/entry)
别名	Apolipoprotein H (YD13515),Apolipoprotein H (YD13515) Rabbit mAb,APOH,APC inhibitor,Activated protein C-binding protein,Anticardiolipin cofactor,Apolipoprotein H,Beta-2-glycoprotein I,B2G1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)