

— ABLYBIO, Help Your Research



# Neutrophil Elastase (YD13499) Rabbit mAb

货号: **AYD15557**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat spleen
细胞定位	Cytoplasmic vesicle, phagosome
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Elastases form a subfamily of serine proteases that hydrolyze many proteins in addition to elastin. Humans have six elastase genes which encode structurally similar proteins. The encoded preproprotein is proteolytically processed to generate the active protease. Following activation, this protease hydrolyzes proteins within specialized neutrophil lysosomes, called azurophil granules, as well as proteins of the extracellular matrix. The enzyme may play a role in degenerative and inflammatory diseases through proteolysis of collagen-IV and elastin. This protein also degrades the outer membrane protein A (OmpA) of <i>E. coli</i> as well as the virulence factors of such bacteria as <i>Shigella</i> , <i>Salmonella</i> and <i>Yersinia</i> . Mutations in this gene are associated with cyclic neutropenia and severe congenital neutropenia (SCN). This gene is present in a gene cluster on chromosome 19.
基因ID	1991
基因名	ELANE
Swiss	P08246 ( <a href="https://www.uniprot.org/uniprotkb/P08246/entry">https://www.uniprot.org/uniprotkb/P08246/entry</a> )
别名	Neutrophil Elastase (YD13499),Neutrophil Elastase (YD13499) Rabbit mAb,ELANE,Bone marrow serine protease,Elastase-2,Human leukocyte elastase,Medullasin,PMN elastase,ELA2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)