

TBX1 (YD13451) Rabbit mAb

货号: **AYD15574**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	Transcription factor that plays a key role in cardiovascular development by promoting pharyngeal arch segmentation during embryonic development (By similarity). Also involved in craniofacial muscle development (By similarity). Together with NKX2-5, acts as a regulator of asymmetric cardiac morphogenesis by promoting expression of PITX2 (By similarity). Acts upstream of TBX1 for the formation of the thymus and parathyroid glands from the third pharyngeal pouch (By similarity). Required for hair follicle stem cell self-renewal (By similarity). Binds to the palindromic T site 5'-TTCACACCTAGGTGTGAA-3' DNA sequence (PubMed:11111039, PubMed:22095455)
基因ID	6899
基因名	TBX1

Swiss	O43435
别名	TBX1 (YD13451)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn