

— ABLYBIO, Help Your Research



BCAR1 (YD13444) Rabbit mAb

货号: **AYD15578**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	93kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cell junction, focal adhesion, Cytoplasm, Cell projection, axon
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	BCAR1, or CAS, is an Src (MIM 190090) family kinase substrate involved in various cellular events, including migration, survival, transformation, and invasion (Sawada et al., 2006 [PubMed 17129785]).[supplied by OMIM, May 2009]
------	---

基因ID	9564
基因名	BCAR1
Swiss	P56945 (https://www.uniprot.org/uniprotkb/P56945/entry)
别名	BCAR1 (YD13444),BCAR1 (YD13444) Rabbit mAb,BCAR1,CRK-associated substrate,Cas scaffolding protein family member 1,p130cas,CAS,CASS1,CRKAS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)