

Kir2.1 (YD13405) Rabbit mAb

货号: **AYD15590**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	48kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse heart,Rat heart
细胞定位	Cell membrane, sarcolemma, T-tubule
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein, which has a greater tendency to allow potassium to flow into a cell rather than out of a cell, probably participates in establishing action potential waveform and excitability of neuronal and muscle tissues. Mutations in this gene have been associated with Andersen syndrome, which is characterized by periodic paralysis, cardiac arrhythmias, and dysmorphic features.
基因ID	3759
基因名	KCNJ2

Swiss	P63252
别名	Kir2.1 (YD13405)

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn