

— ABLYBIO, Help Your Research



# TrkB (YD13377) Rabbit mAb

货号: **AYD15596**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	92kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane, Endosome membrane, Early endosome membrane, Cell projection, axon, dendrite, Cytoplasm, perinuclear region, Postsynaptic density
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the neurotrophic tyrosine receptor kinase (NTRK) family. This kinase is a membrane-bound receptor that, upon neurotrophin binding, phosphorylates itself and members of the M APK pathway. Signalling through this kinase leads to cell differentiation. Mutations in this gene have been associated with obesity and mood disorders. Alternative splicing results in multiple transcript variants.
基因ID	4915
基因名	NTRK2
Swiss	Q16620 ( <a href="https://www.uniprot.org/uniprotkb/Q16620/entry">https://www.uniprot.org/uniprotkb/Q16620/entry</a> )
别名	TrkB (YD13377), TrkB (YD13377) Rabbit mAb, NTRK2, GP145-TrkB, Neurotrophic tyrosine kinase receptor type 2, TrkB tyrosine kinase, Tropomyosin-related kinase B, TRKB

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)