

— ABLYBIO, Help Your Research



PI3-Kinase p85 alpha (YD13226) Rabbit mAb

货号: **AYD15622**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,293F,NIH/3T3
细胞定位	
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Phosphatidylinositol 3-kinase phosphorylates the inositol ring of phosphatidylinositol at the 3-prime position. The enzyme comprises a 110 kD catalytic subunit and a regulatory subunit of either 85, 55, or 50 kD. This gene encodes the 85 kD regulatory subunit. Phosphatidylinositol 3-kinase plays an important role in the metabolic actions of insulin, and a mutation in this gene has been associated with insulin resistance. Alternative splicing of this gene results in four transcript variants encoding different isoforms.
基因ID	5295
基因名	PIK3R1
Swiss	P27986 (https://www.uniprot.org/uniprotkb/P27986/entry)
别名	PI3-Kinase p85 alpha (YD13226),PI3-Kinase p85 alpha (YD13226) Rabbit mAb,PIK3R1,Phosphatidylinositol 3-kinase 85 kDa regulatory subunit alpha,GRB1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)