

— ABLYBIO, Help Your Research



Rab4 (YD13195) Rabbit mAb

货号: **AYD15631**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	24kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Neuro-2a,Mouse brain,Rat brain
细胞定位	Membrane, Cytoplasm, Early endosome membrane, Recycling endosome membrane
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the largest group in the Ras superfamily of small GTPases, which regulate membrane trafficking. The encoded protein is associated with early endosomes and is involved in their sorting and recycling. The protein also plays a role in regulating the recycling of receptors from endosomes to the plasma membrane. Alternatively spliced transcript variants have been observed for this gene.
------	---

基因ID	5867
基因名	RAB4A
Swiss	P20338 (https://www.uniprot.org/uniprotkb/P20338/entry)
别名	Rab4 (YD13195), Rab4 (YD13195) Rabbit mAb, RAB4A, RAB4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)