

— ABLYBIO, Help Your Research



WRN (YD13147) Rabbit mAb

货号: **AYD15645**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	162kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,K-562
细胞定位	Nucleus, nucleolus, nucleoplasm, Chromosome
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the RecQ subfamily and the DEAH (Asp-Glu-Ala-His) subfamily of DNA and RNA helicases. DNA helicases are involved in many aspects of DNA metabolism, including transcription, replication, recombination, and repair. This protein contains a nuclear localization signal in the C-terminus and shows a predominant nucleolar localization. It possesses an intrinsic 3' to 5' DNA helicase activity, and is also a 3' to 5' exonuclease. Based on interactions between this protein and Ku70/80 heterodimer in DNA end processing, this protein may be involved in the repair of double strand DNA breaks. Defects in this gene are the cause of Werner syndrome, an autosomal recessive disorder characterized by premature aging.
基因ID	7486
基因名	WRN
Swiss	Q14191 (https://www.uniprot.org/uniprotkb/Q14191/entry)
别名	WRN (YD13147),WRN (YD13147) Rabbit mAb,WRN,DNA helicase,RecQ-like type 3,RecQ protein-like 2,Werner syndrome protein,DNA 3'-5' helicase WRN,RECQ3,RECQL2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)